



Pulmo-Aide® Compact Compressor Nebulizer System

Quality, Reliability and Affordability

The DeVilbiss Pulmo-Aide Compact Compressor Nebulizer System combines **effective treatment with a small footprint design**. Built to last, the Pulmo-Aide Compact has durable metal mechanical parts and a **five-year warranty** to ensure quality therapy for an extended period of time.

For patients on the go, the Pulmo-Aide Compact offers a convenient carrying handle, **compact design** and **easy transport for AC use**. Reusable nebulizers can be included for convenience, and reliable pressure and flow make it one of the highest quality units on the market.



Pulmo-Aide Compact
Compressor Nebulizer System
(Model 3655D)

Performance and Design Features

- Reliable pressure and flow
- 5-year warranty
- Reusable nebulizer can be included
- Small footprint
- Small, compact design
- Convenient carrying handle
- Metal mechanical parts
- Assembled in USA



Technical Specifications

Built to meet the requirements of HCPCS code:	E0570
Product Weight:	3.8 lbs.
Product Width:	7.5 in.
Product Length/Depth:	7.2 in.
Product Height:	4.0 in.
Shipping Dimensions (box 1):	1 pk. 9.5"w/ 9"d/ 5"h
Shipping Dimensions (box 2):	6 pk. 19.75"w/ 15"d/ 10"h
Ship Method(s):	Ground, Next Day, 2 Day or 3 Day, Freight
Shipping Weight:	1 pk. (5 lbs.) 6 pk. (33 lbs.)
Electrical requirements:	2.5 amps, 115 VAC 60 Hz
Sound level*:	58-62 dBA
Power consumption:	140 watts maximum
Maximum compressor pressure*:	35 psig
Compressor free air flow*:	13 lpm
Operating temperature range:	40 to 104 deg. F (5 to 40 deg. C)
Supplied nebulizer operating pressure*:	12 psig (min)
Supplied nebulizer liter flow*:	5.5 lpm (min)

* Specifications are typical and may vary slightly from system to system.



DeVilbiss Healthcare | 100 DeVilbiss Drive | Somerset, PA 15501 | USA.
800.338.1988 | 814.443.4881 | www.DeVilbissHealthcare.com